

ELECTRIC EXAMINATION TABLE · MOTORIZED HEIGHT · GAS-SPRING BACK PLATE

PCB-MET-200



ABOUT THE PCB-MET-200

Versatile examination, precise positioning.

The **PCB-MET-200 Electric Examination Table** is a versatile examination table for general medical practice, outpatient clinics, gynecology and urology. **Motorized table-top elevation (460–850 mm)** via hand remote provides precise height adjustment between patients.

The **gas-spring back plate** offers strong, smooth recline up to 70° — no motor required, low maintenance. Seat plate, auxiliary plate (with built-in waste basin) and foot pedal are manually adjustable for fast positioning between patients.

The **integrated waterproof mattress** uses foaming technology for comfort and is easy to wipe and disinfect. Available in **Blue, Gray or Black**. A built-in drawer, roll stand and waste basin come standard.

SPECIFICATIONS

Built for daily clinical use.

MODEL PCB-MET-200	LENGTH 1,400 mm	WIDTH 720 mm
TABLETOP LOWEST HEIGHT 460 mm	TABLETOP MAX HEIGHT 850 mm	BACK PLATE UP 70°
SEAT PLATE UP 7°	AUXILIARY PLATE SLIDING 410 mm	FOOT PEDAL ANGLE 15° – 45°
FOOT PEDAL OUTWARD 115°	FOOT PEDAL SLIDING 350 mm	ELEVATION CONTROL Electric motor + hand remote
BACK PLATE MECHANISM Gas-spring	MATTRESS Integrated waterproof	MATTRESS COLORS Blue · Gray · Black
DRAWER Built-in	STANDARD ACCESSORIES Roll stand · Waste basin	OPTIONAL ACCESSORIES Foot switch · Exam light
POWER AC 100 – 240V · 50/60 Hz	CERTIFICATIONS CE · ISO 9001 · ISO 13485	

MULTI-CONTROL ADJUSTMENT

Electric + manual hybrid.

Electric motor handles tabletop elevation while gas-spring and manual controls give clinicians fast, fine adjustment of every other section — combining precision and reliability.



- **Electric:** tabletop height (460–850 mm)
- **Gas spring:** back plate up to 70°
- **Manual:** seat plate, auxiliary plate, foot pedal

- **Hand remote** included (foot switch optional)

COLOR OPTIONS

Choose your mattress color.



Integrated waterproof mattress available in **Blue, Gray** or **Black** — match your clinic décor.

PCB Medical · Kentucky (USA)
Industrial Park Av, Suite 30
Louisville, Kentucky 40041 · Tel. (502) 333-0006

www.pcbmedical.net

Asia Headquarters
PCB Medical · China
2103 Citic Tower · Changsha