

HYDRAULIC OBSTETRIC & GYNECOLOGY TABLE · FOOT-PUMP + HAND-WHEEL · STAINLESS STEEL · NO POWER

PCB-OMT-250



ABOUT THE PCB-OMT-250

Reliable hydraulic positioning, no electricity required.

The **PCB-OMT-250** is a fully hydraulic **obstetric and gynecology table**. **Hand-wheel and foot-pump operation** provide precise positioning without electrical dependence — ideal for facilities with unstable power. The **polished stainless steel frame** resists corrosion and stands up to daily cleaning protocols.

MANUAL HYDRAULIC

Manual hydraulic lift & positioning.



Foot-pump elevation plus **hand-wheel positioning** gives clinicians fine control of every section — back, leg, and tabletop tilt — without depending on the building's power supply.

- Height range **750 – 980 mm** via foot pump
- Back section **up to 80°** via hand wheel
- Leg section **down to 90°** below horizontal
- Trendelenburg / reverse Trendelenburg **±15°**
- Load capacity **200 kg**

SPECIFICATIONS

Complete technical data.

MODEL

PCB-OMT-250

TABLETOP LENGTH

1,850 mm

TABLETOP WIDTH

550 mm

HEIGHT RANGE

750 – 980 mm

BACK SECTION UP

0° – 80°

LEG SECTION DOWN

0° – 90°

TRENDELENBURG / REVERSE

±15°

LOAD CAPACITY

200 kg

LIFT MECHANISM

Hydraulic (manual)

FRAME

Stainless steel

UPHOLSTERY

Seamless PU leather

CASTORS

4 x Ø125 mm with brake

REMOVABLE BASIN

Included

ACCESSORIES

Stirrups + leg supports

CERTIFICATIONS

CE · ISO 9001 · ISO 13485

PCB Medical · Kentucky (USA)

Industrial Park Av, Suite 30
Louisville, Kentucky 40041 · Tel. (502) 333-0006

www.pcbmedical.net

Asia Headquarters

PCB Medical · China
2103 Citic Tower · Changsha