



SKLAR OR QUALITY · PRE-CONFIGURED SET

AMPUTATION SET

Pre-configured 73-piece instrument set by Sklar Corporation (ref. 98-1017) — German stainless steel, factory-inspected, backed by Sklar's Lifetime Guarantee against manufacturing defects. Distributed in your market by PCB Medical.

German Stainless Steel

Lifetime Guarantee

51-point QC

5-yr carbide warranty

Autoclave-ready

Set composition

Standard PCB Medical / Sklar 98-1017 configuration. Custom configurations available on request.

QTY	INSTRUMENT	SKLAR REF.
2	Sklar® Scalpel Handle - #3	06-2903
1	Sklar® Scalpel Handle - #4	06-2904
1	Sklar Edge™ TC Mayo Dissecting Scissors - 6-3/4", Straight	16-1510
1	Sklar Edge™ TC Mayo Dissecting Scissors - 6-3/4", Curved	16-1610
1	Sklar Edge™ TC Metzenbaum Dissecting Scissors - 7", Curved, Blunt/Blunt	16-2110
6	Sklar® Rochester-Pean Hemostatic Forceps - 6-1/4", Curved	17-2262
2	Sklar® Rochester-Pean Hemostatic Forceps - 8", Curved	17-2280
4	Sklar® Rochester-Ochsner Forceps - 6-1/4", Straight, Serrated, 1×2 Teeth	17-2862
2	Sklar® Rochester-Ochsner Forceps - 8", Straight, Serrated, 1×2 Teeth	17-2880
12	Sklar® Crile Hemostatic Forceps - 5-1/2", Curved	17-3155
2	Sklar® Dressing Forceps - 6", Straight, Serrated	19-1060
1	Sklar® Tissue Forceps - 6", Straight, 1×2 Teeth	19-1260
1	Sklar® TC Halsey / Webster Needle Holder - 5", Serrated	21-8001
2	Sklar® TC Mayo-Hegar Needle Holder - 6", Serrated	21-8031
2	Sklar® TC Mayo-Hegar Needle Holder - 7", Serrated	21-8032
4	Sklar® Allis Tissue Forceps - 6", 5×6 Teeth"	36-2260
1	Sklar® Liston Bone Cutting Forceps	40-3055
1	Sklar® Liston Bone Cutting Forceps	40-3077
1	Sklar® Stille-Liston Bone Cutting Forceps	40-3810

73 Total pieces in this set (37 distinct instruments)

QTY	INSTRUMENT	SKLAR REF.
1	Sklar® Stille-Luer Rongeur	40-4085
1	Sklar® Leksell Rongeur	40-4092
1	Sklar® Leksell Rongeur	40-4095
1	Sklar® Beyer Rongeur - 7", Curved, 3mm, Double-Action	40-4268
1	Sklar® Solid Bar Handle For Gigli Saw	40-5520
1	Sklar® Percy Amputation Retractor	40-5551
1	Sklar® Bone Hook	40-6857
1	Sklar® Key Elevator	40-7402
1	Sklar® Key Elevator	40-7404
1	Sklar® Key Elevator	40-7406
1	Sklar® Putti Bone Raspatory	40-9812
2	Sklar® Joseph Skin Hook - 6-1/4", 1 Prong, Sharp	47-2500
4	Sklar® Backhaus Towel Clamp - 5-1/4"	47-2955
2	Sklar® Adson Tissue Forceps - 4-3/4", Straight, 1×2 Teeth	50-3047
1	Sklar® US Army Retractor	60-1085
2	Sklar® Volkman Retractor - 8-1/2", Sharp, 6 Prongs, Ring Handle	60-5826
2	Sklar® Senn Retractor - 6-1/4", Sharp Prongs	60-6162
2	Sklar® Foerster Sponge Forceps - 9-1/2", Straight, Serrated	87-2195
73	Total pieces in this set (37 distinct instruments)	



Sklar — A 130-year heritage of surgical excellence

Sklar Corporation has been crafting surgical instruments in West Chester, Pennsylvania since 1882. Every instrument is forged from high-grade German stainless steel, finished by hand, and verified by Sklar's documented quality system before shipment.



51-point incoming inspection

Each component is checked against Sklar's quality standards before being approved for assembly.



6-point final inspection

Every assembled instrument is hand-checked against a final 6-point checklist before being shipped.

DE

German stainless steel

Surgical-grade material on OR Quality and Sklarlite™ lines. Hand-finished to exacting tolerances.



130+ years heritage

Crafting surgical instruments since 1882. Trusted by U.S. hospitals, surgery centers, and clinics.

Set reference data

SKLAR REFERENCE	98-1017
SET TYPE	Amputation Set — pre-configured (73 pieces)
QUALITY TIER	Sklar OR Quality — German stainless
REUSABILITY	Reusable · non-sterile · steam autoclave compatible
INSPECTION	51-point incoming + 6-point final (per Sklar IFU)
WARRANTY · STAINLESS	Lifetime against manufacturing defects
WARRANTY · CARBIDE	5 years (sharpen / repair / replace)
MANUFACTURER	Sklar Corporation · West Chester, PA · est. 1882
DISTRIBUTION	PCB Medical — authorized Sklar distributor

WARRANTY & QUALITY

Backed by Sklar's Lifetime Guarantee

"All Sklar Instruments are guaranteed for life against manufacturing defects, provided that the instrument is used for its intended surgical purpose. Sklar will repair or replace defective instruments at no charge, at our discretion."

✓ Lifetime guarantee

✓ 5-yr carbide warranty

✓ 51-point QC

✓ Made by Sklar

