



SKLAR OR QUALITY · PRE-CONFIGURED SET

LAMINECTOMY SET

Pre-configured 112-piece instrument set by Sklar Corporation (ref. 98-1032) — German stainless steel, factory-inspected, backed by Sklar's Lifetime Guarantee against manufacturing defects. Distributed in your market by PCB Medical.

German Stainless Steel

Lifetime Guarantee

51-point QC

5-yr carbide warranty

Autoclave-ready

Set composition

Standard PCB Medical / Sklar 98-1032 configuration. Custom configurations available on request.

QTY	INSTRUMENT	SKLAR REF.
1	Sklar® Scalpel Handle - #3	06-2903
1	Sklar® Scalpel Handle - #4	06-2904
1	Sklar® Scalpel Handle - #7	06-2907
1	Sklar® Operating Scissors - 5-1/2", Straight, Sharp/Blunt	14-1055
1	Sklar® Operating Scissors - 5-1/2", Curved, Sharp/Blunt	14-2055
1	Sklar® Mayo Dissecting Scissors - 6-3/4", Straight	15-1567
1	Sklar® Mayo Dissecting Scissors - Curved, 6-3/4"	15-2567
12	Sklar® Halsted Mosquito Hemostatic Forceps - 5", Straight	17-1450
12	Sklar® Halsted Mosquito Hemostatic Forceps - 5", Curved	17-1550
12	Sklar® Crile Forceps - 6-1/4"	17-3062
12	Sklar® Crile Forceps - 6-1/4"	17-3162
1	Sklar® Dressing Forceps - 5-1/2", Straight, Serrated	19-1055
2	Sklar® Dressing Forceps - 8", Straight, Serrated	19-1080
1	Sklar® Tissue Forceps - 5-1/2", Straight, 1×2 Teeth	19-1255
2	Sklar® Tissue Forceps - 8", Straight, 1×2 Teeth	19-1280
1	Sklar® Cushing Brain Forceps	22-6770
1	Sklar® Cushing Brain Forceps	22-6870
4	Sklar® Allis Tissue Forceps - 6", 5×6 Teeth"	36-2260
1	Sklar® Luer Bone Rongeur	40-3670

112 Total pieces in this set (45 distinct instruments)

QTY	INSTRUMENT	SKLAR REF.
1	Sklar® Stille-Liston Bone Cutting Forceps	40-3710
1	Sklar® Stille-Horsley Bone Cutting Forceps	40-3910
1	Sklar® Stille-Luer Rongeur	40-3985
1	Sklar® Stille-Luer Rongeur	40-4090
1	Sklar® Leksell Rongeur	40-4092
1	Sklar® Bruns Bone Curette - 6-3/4", Straight, Oval, #000	40-7000
1	Sklar® Bruns Bone Curette - 6-3/4", Straight, Oval, #00	40-7100
1	Sklar® Bruns Bone Curette - 6-3/4", Straight, Oval, #0	40-7110
1	Sklar® Bruns Bone Curette - 6-3/4", Straight, Oval, #1	40-7111
1	Sklar® Bruns Bone Curette - 6-3/4", Straight, Oval, #2	40-7112
1	Sklar® Bruns Bone Curette - 6-3/4", Straight, Oval, #3	40-7113
1	Sklar® Bruns Bone Curette - 6-3/4", Straight, Oval, #4	40-7114
1	Sklar® Bruns Bone Curette - 6-3/4", Straight, Oval, #5	40-7115
1	Sklar® Bruns Bone Curette - 6-3/4", Straight, Oval, #6	40-7116
2	Sklar® Weitlaner Retractor - Sharp 3×4 Teeth 8"	42-1975
1	Sklar® Love-Gruenwald Intervertebral Disc Rongeur	42-3150
1	Sklar® Love-Gruenwald Intervertebral Disc Rongeur	42-3156
1	Sklar® Love-Gruenwald Intervertebral Disc Rongeur	42-3158
1	Sklar® Love-Gruenwald Intervertebral Disc Rongeur	42-3160
1	Sklar® Love-Kerrison Laminectomy Rongeur	42-3418
1	Sklar® Love-Kerrison Laminectomy Rongeur	42-3429
12	Sklar® Backhaus Towel Clamp - 5-1/4"	47-2955
1	Sklar® Adson Bayonet Dressing Forceps - 7-1/2", Straight, Serrated, Delicate	50-3371
2	Sklar® Beckman-Adson Retractor - 12-1/2", Sharp	57-1322
1	Sklar® Metzenbaum Dissecting Scissors - 7", Straight, Blunt/Blunt	75-5570
6	Sklar® Foerster Sponge Forceps - 9-1/2", Straight, Serrated	87-2195

112 Total pieces in this set (45 distinct instruments)



Sklar — A 130-year heritage of surgical excellence

Sklar Corporation has been crafting surgical instruments in West Chester, Pennsylvania since 1882. Every instrument is forged from high-grade German stainless steel, finished by hand, and verified by Sklar's documented quality system before shipment.



51-point incoming inspection

Each component is checked against Sklar's quality standards before being approved for assembly.



6-point final inspection

Every assembled instrument is hand-checked against a final 6-point checklist before being shipped.

DE

German stainless steel

Surgical-grade material on OR Quality and Sklarlite™ lines. Hand-finished to exacting tolerances.



130+ years heritage

Crafting surgical instruments since 1882. Trusted by U.S. hospitals, surgery centers, and clinics.

Set reference data

SKLAR REFERENCE	98-1032
SET TYPE	Laminectomy Set — pre-configured (112 pieces)
QUALITY TIER	Sklar OR Quality — German stainless
REUSABILITY	Reusable · non-sterile · steam autoclave compatible
INSPECTION	51-point incoming + 6-point final (per Sklar IFU)
WARRANTY · STAINLESS	Lifetime against manufacturing defects
WARRANTY · CARBIDE	5 years (sharpen / repair / replace)
MANUFACTURER	Sklar Corporation · West Chester, PA · est. 1882
DISTRIBUTION	PCB Medical — authorized Sklar distributor

WARRANTY & QUALITY

Backed by Sklar's Lifetime Guarantee

"All Sklar Instruments are guaranteed for life against manufacturing defects, provided that the instrument is used for its intended surgical purpose. Sklar will repair or replace defective instruments at no charge, at our discretion."

✓ Lifetime guarantee

✓ 5-yr carbide warranty

✓ 51-point QC

✓ Made by Sklar

