



SKLAR OR QUALITY · PRE-CONFIGURED SET

# PACEMAKER SET

Pre-configured 72-piece instrument set by Sklar Corporation (ref. 98-1043) — German stainless steel, factory-inspected, backed by Sklar's Lifetime Guarantee against manufacturing defects. Distributed in your market by PCB Medical.

German Stainless Steel

Lifetime Guarantee

51-point QC

5-yr carbide warranty

Autoclave-ready

## Set composition

Standard PCB Medical / Sklar 98-1043 configuration. Custom configurations available on request.

QTY	INSTRUMENT	SKLAR REF.
2	Sklar® Scalpel Handle - #3	06-2903
1	Sklar® Scalpel Handle - #7	06-2907
1	Sklar Edge™ TC Mayo Dissecting Scissors - 6-3/4", Straight	16-1510
1	Sklar Edge™ TC Mayo Dissecting Scissors - 6-3/4", Curved	16-1610
1	Sklar Edge™ TC Metzenbaum-Lahey Dissecting Scissors - 5-3/4", Curved, Blunt/Blunt, Delicate	16-1905
1	Sklar Edge™ TC Metzenbaum Dissecting Scissors - 7", Curved, Blunt/Blunt	16-2110
4	Sklar® Halsted Mosquito Hemostatic Forceps - 5", Curved	17-1550
2	Sklar® Mixer Forceps Right Angle - 6-1/4"	17-1762
4	Sklar® Rochester-Pean Hemostatic Forceps - 6-1/4", Curved	17-2262
4	Sklar® Rochester-Ochsner Forceps - 6-1/4", Straight, Serrated, 1×2 Teeth	17-2862
4	Sklar® Crile Hemostatic Forceps - 5-1/2", Straight	17-3055
10	Sklar® Crile Hemostatic Forceps - 5-1/2", Curved	17-3155
1	Sklar® Dressing Forceps - 6", Straight, Serrated	19-1060
2	Sklar® Tissue Forceps - 6", Straight, 1×2 Teeth	19-1260
2	Sklar® TC Halsey / Webster Needle Holder - 5", Serrated	21-8001
2	Sklar® TC Mayo-Hegar Needle Holder - 6", Serrated	21-8031
1	Sklar® Weitlaner Retractor - Sharp 3×4 Teeth 5-1/2"	22-8255
2	Sklarlite XD™ Mixer Right Angle Forceps	23-2521
<b>72</b>	<b>Total pieces in this set (30 distinct instruments)</b>	

QTY	INSTRUMENT	SKLAR REF.
4	Sklar® Allis Tissue Forceps - 6", 4×5 Teeth	36-2160
2	Sklar® Backhaus Towel Clamp - 5-1/4"	47-2955
2	Sklar® Cushing Nerve and Vein Retractor	50-1290
4	Sklar® Adson Hemostatic Forceps - 7-1/4", Curved, Half-Serrated Jaws	50-1472
2	Sklar® Adson Tissue Forceps - 4-3/4", Straight, 1×2 Teeth	50-3047
2	Sklar® Adson (Ewald) Dressing Forceps - 4-3/4", Straight, Serrated	50-3147
2	Sklar® DeBakey Atraumatic Tissue Forceps	52-5162
1	Sklar® US Army Retractor	60-1085
2	Sklar® Kelly Retractor	60-1656
2	Sklar® Richardson Retractor	60-1676
2	Sklar® Volkman Retractor - 8-1/2", Sharp, 4 Prongs, Ring Handle	60-5824
2	Sklar® Foerster Sponge Forceps - 9-1/2", Straight, Serrated	87-2195
<b>72</b>	<b>Total pieces in this set (30 distinct instruments)</b>	



## Sklar — A 130-year heritage of surgical excellence

Sklar Corporation has been crafting surgical instruments in West Chester, Pennsylvania since 1882. Every instrument is forged from high-grade German stainless steel, finished by hand, and verified by Sklar's documented quality system before shipment.



### 51-point incoming inspection

Each component is checked against Sklar's quality standards before being approved for assembly.



### 6-point final inspection

Every assembled instrument is hand-checked against a final 6-point checklist before being shipped.

### DE

#### German stainless steel

Surgical-grade material on OR Quality and Sklarlite™ lines. Hand-finished to exacting tolerances.



### 130+ years heritage

Crafting surgical instruments since 1882. Trusted by U.S. hospitals, surgery centers, and clinics.

## Set reference data

SKLAR REFERENCE	98-1043
SET TYPE	Pacemaker Set — pre-configured (72 pieces)
QUALITY TIER	Sklar OR Quality — German stainless
REUSABILITY	Reusable · non-sterile · steam autoclave compatible
INSPECTION	51-point incoming + 6-point final (per Sklar IFU)
WARRANTY · STAINLESS	Lifetime against manufacturing defects
WARRANTY · CARBIDE	5 years (sharpen / repair / replace)
MANUFACTURER	Sklar Corporation · West Chester, PA · est. 1882
DISTRIBUTION	PCB Medical — authorized Sklar distributor

### WARRANTY & QUALITY

## Backed by Sklar's Lifetime Guarantee

*"All Sklar Instruments are guaranteed for life against manufacturing defects, provided that the instrument is used for its intended surgical purpose. Sklar will repair or replace defective instruments at no charge, at our discretion."*

✓ Lifetime guarantee

✓ 5-yr carbide warranty

✓ 51-point QC

✓ Made by Sklar

