



SKLAR OR QUALITY · PRE-CONFIGURED SET

BREAST AUGMENT/REDUCTION SET

Pre-configured 54-piece instrument set by Sklar Corporation (ref. 98-1095) — German stainless steel, factory-inspected, backed by Sklar's Lifetime Guarantee against manufacturing defects. Distributed in your market by PCB Medical.

German Stainless Steel

Lifetime Guarantee

51-point QC

5-yr carbide warranty

Autoclave-ready

Set composition

Standard PCB Medical / Sklar 98-1095 configuration. Custom configurations available on request.

QTY	INSTRUMENT	SKLAR REF.
1	Steel Ruler - 6"	06-2860
2	Sklar® Scalpel Handle - #3R, Straight, Round	06-2908
1	Emesis Basin - 6-5/8" x 3-1/2" x 1-1/2", 10 fl oz	10-1497
1	Sklar® Mayo Dissecting Scissors - Curved, 6-3/4"	15-2567
4	Sklar® Hartmann Mosquito Forceps - 3-1/2", Straight	17-1235
2	Sklar® Halsted Mosquito Hemostatic Forceps - 5", Straight	17-1450
2	Sklar® Rochester-Ochsner (Kocher) Forceps - 5-1/2", Straight, Serrated, 1×2 Teeth	17-2855
2	Sklar® Tissue Forceps - 4-1/2", Straight, 1×2 Teeth	19-1245
2	Sklar® Webster / Halsey Needle Holder - 5-1/4", Smooth	20-1350
2	Sklar® Baumgartner Needle Holder - 5-1/4", Serrated	20-1450
3	Sklar® Olsen-Hegar Needle Holder - 5-1/2", Serrated	20-2755
1	Sklar® Heaney Needle Holder - 8-1/4", Curved, Serrated	20-4682
1	Sklar® Goldman-Fox Scissors - 5", Curved, Serrated, Sharp/Sharp	22-2520
2	Sklar® Allis Tissue Forceps - 6", 4×5 Teeth	36-2160
1	Sklar® Iris Scissors - 4-1/2", Straight, Sharp/Sharp	47-1145
1	Sklar® Iris Scissors - 4-1/2", Curved, Sharp/Sharp	47-1245
1	Sklar® Joseph Skin Hook - 6-1/4", 2 Prongs, Sharp, 2mm	47-2502
1	Sklar® Joseph Skin Hook - 6-1/4", 2 Prongs, Sharp, 5mm	47-2505
4	Sklar® Backhaus Towel Clamp - 3-1/2"	47-2935

54 Total pieces in this set (34 distinct instruments)

QTY	INSTRUMENT	SKLAR REF.
1	Sklar® Frazier Suction Tube - 8 French, Angled 30°	50-2008
2	Sklar® Adson Tissue Forceps - 4-3/4", Straight, 1×2 Teeth	50-3047
1	Sklar® Adson (Ewald) Dressing Forceps - 4-3/4", Straight, Serrated	50-3147
1	Sklar® DeBakey Atraumatic Tissue Forceps	52-5177
2	Sklar® US Army Retractor	60-1085
2	Sklar® Deaver #4 2X12	60-3312
3	Sklar® Zalkind Ribbon Retractor - 7-3/4" x 1/2"	60-3596
1	Sklar® Senn Retractor - 6-1/4", Sharp Prongs	60-6162
1	Sklar® Stevens Tenotomy Scissors - 4-1/2", Straight, Blunt/Blunt	64-1242
1	Sklar® Stevens Tenotomy Scissors - 4-1/2", Curved, Blunt/Blunt	64-1342
1	Sklar® Bowman Lacrimal Probe	65-4107
1	Sklar® Schnidt Tonsil Forceps	74-3075
1	Sklar® Metzenbaum Dissecting Scissors - 7", Straight, Blunt/Blunt	75-5570
1	Sklar® Metzenbaum Dissecting Scissors - 7", Curved, Blunt/Blunt	75-5670
1	Sklar® Electrosurgical Adson Bipolar Tissue Forceps	91-6100
54	Total pieces in this set (34 distinct instruments)	



Sklar — A 130-year heritage of surgical excellence

Sklar Corporation has been crafting surgical instruments in West Chester, Pennsylvania since 1882. Every instrument is forged from high-grade German stainless steel, finished by hand, and verified by Sklar's documented quality system before shipment.



51-point incoming inspection

Each component is checked against Sklar's quality standards before being approved for assembly.



6-point final inspection

Every assembled instrument is hand-checked against a final 6-point checklist before being shipped.

DE

German stainless steel

Surgical-grade material on OR Quality and Sklarlite™ lines. Hand-finished to exacting tolerances.



130+ years heritage

Crafting surgical instruments since 1882. Trusted by U.S. hospitals, surgery centers, and clinics.

Set reference data

SKLAR REFERENCE	98-1095
SET TYPE	Breast Augment/Reduction Set — pre-configured (54 pieces)
QUALITY TIER	Sklar OR Quality — German stainless
REUSABILITY	Reusable · non-sterile · steam autoclave compatible
INSPECTION	51-point incoming + 6-point final (per Sklar IFU)
WARRANTY · STAINLESS	Lifetime against manufacturing defects
WARRANTY · CARBIDE	5 years (sharpen / repair / replace)
MANUFACTURER	Sklar Corporation · West Chester, PA · est. 1882
DISTRIBUTION	PCB Medical — authorized Sklar distributor

WARRANTY & QUALITY

Backed by Sklar's Lifetime Guarantee

"All Sklar Instruments are guaranteed for life against manufacturing defects, provided that the instrument is used for its intended surgical purpose. Sklar will repair or replace defective instruments at no charge, at our discretion."

✓ Lifetime guarantee

✓ 5-yr carbide warranty

✓ 51-point QC

✓ Made by Sklar

