



SKLAR OR QUALITY · PRE-CONFIGURED SET

RHYTIDECTOMY/FACELIFT SET

Pre-configured 42-piece instrument set by Sklar Corporation (ref. 98-1097) — German stainless steel, factory-inspected, backed by Sklar's Lifetime Guarantee against manufacturing defects. Distributed in your market by PCB Medical.

German Stainless Steel

Lifetime Guarantee

51-point QC

5-yr carbide warranty

Autoclave-ready

Set composition

Standard PCB Medical / Sklar 98-1097 configuration. Custom configurations available on request.

QTY	INSTRUMENT	SKLAR REF.
1	Steel Ruler - 6"	06-2860
1	Sklar® Scalpel Handle - #7	06-2907
1	Sklar® Scalpel Handle - #3, with Metric Graduations	06-2910
1	Skларcut™ Castanares Face Lift Scissors	15-3549
1	Skларcut™ Iris Scissors - 4-1/2", Straight, Sharp/Sharp, 1 Standard Edge, 1 Micro-Serrated Edge	15-3585
1	Skларcut™ Iris Scissors - 4-1/2", Curved, Sharp/Sharp, 1 Standard Edge, 1 Micro-Serrated Edge	15-3590
1	Skлар Edge™ TC Mayo Dissecting Scissors - 6-3/4", Straight	16-1510
1	Skлар Edge™ TC Metzenbaum-Lahey Dissecting Scissors - 5-3/4", Curved, Blunt/Blunt, Delicate	16-1905
2	Skлар® Halsted Mosquito Hemostatic Forceps - 5", Straight	17-1450
2	Skлар® Halsted Mosquito Hemostatic Forceps - 5", Curved	17-1550
1	Skлар® TC Webster / Halsey Needle Holder - 5-1/4", Smooth	21-8002
1	Skлар® TC Crile-Wood / DeBakey Needle Holder - 6", Serrated, Narrow Jaws	21-8006
2	Skлар® TC Webster / Halsey Needle Holder - 5-1/4", Serrated	21-8080
1	Skлар® Tissue Forceps - 5-1/2", Straight, 1x2 Teeth, Delicate	22-1063
1	Skлар® Dressing Forceps - 5", Straight, Serrated, Delicate	22-1162
1	Skлар® Ragnell Scissors - 5" Straight	22-2653
1	Skлар® Spencer (Littauer) Suture Scissors - 3-1/2"	22-2835

42 Total pieces in this set (33 distinct instruments)

QTY	INSTRUMENT	SKLAR REF.
2	Sklar® Allis Tissue Forceps - 6", 4×5 Teeth	36-2160
1	Sklar® Rhytidectomy Scissors	41-1105
1	Sklar® Blepharoplasty Scissors - 4-3/4"	47-1147
1	Sklar® Adson-Brown Tissue Forceps - 4-3/4", Straight, Side-Grasping 7×7 Teeth	47-1447
2	Sklar® Joseph Skin Hook - 6-1/4", 2 Prongs, Sharp, 10mm	47-2510
2	Sklar® Backhaus Towel Clamp - 3-1/2"	47-2935
2	Sklar® Backhaus Towel Clamp - 5-1/4"	47-2955
1	Sklar® Frazier Suction Tube - 8 French, Angled 30°	50-2008
1	Sklar® Adson Tissue Forceps - 4-3/4", Straight, 1×2 Teeth	50-3047
1	Sklar® Potts-Smith Tissue Forceps - 7", Straight, 1×2 Teeth	52-3370
1	Sklar® DeBakey Atraumatic Forceps - 6-1/4"	52-5062
1	Sklar® DeBakey Atraumatic Tissue Forceps	52-5162
2	Sklar® US Army Retractor	60-1085
1	Sklar® Richardson Retractor	60-1676
1	Sklar® DEAVER SKL GR HDL #3 1-1/2 X12	60-2812
2	Sklar® Senn Retractor - 6-1/4", Sharp Prongs	60-6162
42	Total pieces in this set (33 distinct instruments)	



Sklar — A 130-year heritage of surgical excellence

Sklar Corporation has been crafting surgical instruments in West Chester, Pennsylvania since 1882. Every instrument is forged from high-grade German stainless steel, finished by hand, and verified by Sklar's documented quality system before shipment.



51-point incoming inspection

Each component is checked against Sklar's quality standards before being approved for assembly.



6-point final inspection

Every assembled instrument is hand-checked against a final 6-point checklist before being shipped.

DE

German stainless steel

Surgical-grade material on OR Quality and Sklarlite™ lines. Hand-finished to exacting tolerances.



130+ years heritage

Crafting surgical instruments since 1882. Trusted by U.S. hospitals, surgery centers, and clinics.

Set reference data

SKLAR REFERENCE	98-1097
SET TYPE	Rhytidectomy/Facelift Set — pre-configured (42 pieces)
QUALITY TIER	Sklar OR Quality — German stainless
REUSABILITY	Reusable · non-sterile · steam autoclave compatible
INSPECTION	51-point incoming + 6-point final (per Sklar IFU)
WARRANTY · STAINLESS	Lifetime against manufacturing defects
WARRANTY · CARBIDE	5 years (sharpen / repair / replace)
MANUFACTURER	Sklar Corporation · West Chester, PA · est. 1882
DISTRIBUTION	PCB Medical — authorized Sklar distributor

WARRANTY & QUALITY

Backed by Sklar's Lifetime Guarantee

"All Sklar Instruments are guaranteed for life against manufacturing defects, provided that the instrument is used for its intended surgical purpose. Sklar will repair or replace defective instruments at no charge, at our discretion."

✓ Lifetime guarantee

✓ 5-yr carbide warranty

✓ 51-point QC

✓ Made by Sklar

