

SKLAR OR QUALITY · PRE-CONFIGURED SET

FST MINOR BASIC INSTRU SET

Pre-configured 67-piece instrument set by Sklar Corporation (ref. 98-1324) — German stainless steel, factory-inspected, backed by Sklar's Lifetime Guarantee against manufacturing defects. Distributed in your market by PCB Medical.

German Stainless Steel

Lifetime Guarantee

51-point QC

5-yr carbide warranty

Autoclave-ready

Set composition

Standard PCB Medical / Sklar 98-1324 configuration. Custom configurations available on request.

QTY	INSTRUMENT	SKLAR REF.
2	Sklar® Frazier Suction Tube - 8 French, Angled 30°	50-2008
2	Sklar® Yankauer Suction Tube	07-1801
2	Sklar® Gerald Dressing Forceps - 7", Straight, Serrated	47-1180
2	Sklar® Foerster Sponge Forceps - 9-1/2", Straight, Serrated	87-2195
6	Sklar® Crile Hemostatic Forceps - 5-1/2", Curved	17-3155
6	Sklar® Halsted Mosquito Hemostatic Forceps - 5", Curved	17-1550
2	Sklar® Mixer Forceps Right Angle - 7-1/4"	17-1772
4	Sklar® Rochester-Ochsner Forceps - 6-1/4", Straight, Serrated, 1x2 Teeth	17-2862
4	Sklar® Rochester-Pean Hemostatic Forceps - 6-1/4", Curved	17-2262
2	Sklar® Tissue Forceps - 5-1/2", Straight, 1x2 Teeth	19-1255
2	Sklar® Adson Tissue Forceps - 4-3/4", Straight, 1x2 Teeth, with Smooth Tying Platform	50-3049
4	Sklar® Allis Tissue Forceps - 6", 4x5 Teeth	36-2160
2	Sklar® DeBakey Atraumatic Tissue Forceps	52-5177
6	Sklar® Backhaus Towel Clamp - 5-1/4"	47-2955
2	Sklar® Scalpel Handle - #7	06-2907
1	Sklar® TC DeBakey Cardiovascular Needle Holder - 7", Serrated	21-8087
2	Sklar® TC Mayo-Hegar Needle Holder - 7", Serrated	21-8032
1	Sklar® TC Ryder Needle Holder - 5", Serrated	21-8019
1	Sklar® Weinstein Instrument Holder	10-2012

67 Total pieces in this set (31 distinct instruments)

QTY	INSTRUMENT	SKLAR REF.
2	Sklar® Weitlaner Retractor - Blunt 3x4 Teeth 6-1/2"	42-1865
2	Sklar® Gelpi Retractor	22-8855
1	Sklar® US Army Retractor	60-1085
1	Sklar® Hupp Tracheal Retractor	78-1113
1	Sklar® Cushing Nerve and Vein Retractor	50-1290
1	Sklar® Mayo Dissecting Scissors - Curved, 6-3/4"	15-2567
1	Sklar® Mayo Dissecting Scissors - 6-3/4", Straight	15-1567
1	Sklar® Metzenbaum Dissecting Scissors - 7", Curved, Blunt/Blunt	75-5670
1	Sklar® Metzenbaum-Nelson Bariatric Dissecting Scissors - 11", Curved, Blunt/Blunt	55-9411
1	Sklar® Iris Scissors - 4-1/2", Straight, Sharp/Sharp	47-1145
1	Sklar® Metzenbaum-Lahey Dissecting Scissors - 5-3/4", Curved, Blunt/Blunt	47-1360
1	Sklar® Stevens Tenotomy Scissors - 4-1/2", Curved, Blunt/Blunt	64-1342
67	Total pieces in this set (31 distinct instruments)	



Sklar — A 130-year heritage of surgical excellence

Sklar Corporation has been crafting surgical instruments in West Chester, Pennsylvania since 1882. Every instrument is forged from high-grade German stainless steel, finished by hand, and verified by Sklar's documented quality system before shipment.



51-point incoming inspection

Each component is checked against Sklar's quality standards before being approved for assembly.



6-point final inspection

Every assembled instrument is hand-checked against a final 6-point checklist before being shipped.

DE

German stainless steel

Surgical-grade material on OR Quality and Sklarlite™ lines. Hand-finished to exacting tolerances.



130+ years heritage

Crafting surgical instruments since 1882. Trusted by U.S. hospitals, surgery centers, and clinics.

Set reference data

SKLAR REFERENCE	98-1324
SET TYPE	Fst Minor Basic Instru Set — pre-configured (67 pieces)
QUALITY TIER	Sklar OR Quality — German stainless
REUSABILITY	Reusable · non-sterile · steam autoclave compatible
INSPECTION	51-point incoming + 6-point final (per Sklar IFU)
WARRANTY · STAINLESS	Lifetime against manufacturing defects
WARRANTY · CARBIDE	5 years (sharpen / repair / replace)
MANUFACTURER	Sklar Corporation · West Chester, PA · est. 1882
DISTRIBUTION	PCB Medical — authorized Sklar distributor

WARRANTY & QUALITY

Backed by Sklar's Lifetime Guarantee

"All Sklar Instruments are guaranteed for life against manufacturing defects, provided that the instrument is used for its intended surgical purpose. Sklar will repair or replace defective instruments at no charge, at our discretion."

✓ Lifetime guarantee

✓ 5-yr carbide warranty

✓ 51-point QC

✓ Made by Sklar

