

SKLAR OR QUALITY · PRE-CONFIGURED SET

MINOR BASIC SET

Pre-configured 41-piece instrument set by Sklar Corporation (ref. 98-1331) — German stainless steel, factory-inspected, backed by Sklar's Lifetime Guarantee against manufacturing defects. Distributed in your market by PCB Medical.

German Stainless Steel

Lifetime Guarantee

51-point QC

5-yr carbide warranty

Autoclave-ready

Set composition

Standard PCB Medical / Sklar 98-1331 configuration. Custom configurations available on request.

QTY	INSTRUMENT	SKLAR REF.
4	Sklar® Backhaus Towel Clamp - 3-1/2"	47-2935
2	Sklar® Halsted Mosquito Hemostatic Forceps - 5", Curved	17-1550
4	Sklar® Kelly Hemostatic Forceps - 5-1/2", Curved	17-2155
2	Sklar® Allis Tissue Forceps - 6", 3×4 Teeth	36-2060
2	Sklar® Rochester-Ochsner (Kocher) Forceps - 5-1/2", Straight, Serrated, 1×2 Teeth	17-2855
2	Sklar® Mixer Forceps	22-3270
2	Sklar® Ballenger Sponge Forceps - 7", Straight, Smooth	74-1370
1	Sklar® TC Crile-Wood / DeBakey Needle Holder - 6", Serrated, Narrow Jaws	21-8006
1	Sklar® TC Webster / Halsey Needle Holder - 5-1/4", Smooth	21-8002
1	Sklar® Mayo Dissecting Scissors - 6-3/4", Straight	15-1567
1	Sklar® Metzenbaum-Lahey Dissecting Scissors - 5-3/4", Curved, Blunt/Blunt	47-1360
1	Sklar® Metzenbaum-Lahey Dissecting Scissors - 4-1/2", Curved, Blunt/Blunt	22-2545
2	Sklar® Scalpel Handle - #3, with Metric Graduations	06-2910
1	Sklar® DeBakey Atraumatic Tissue Forceps	52-5162
1	Sklar® DeBakey Atraumatic Tissue Forceps	52-5177
2	Sklar® Adson Tissue Forceps - 4-3/4", Straight, 1×2 Teeth, with Smooth Tying Platform	50-3049
2	Sklar® Adson-Brown Tissue Forceps - 4-3/4", Straight, Side-Grasping 7×7 Teeth	47-1447
2	Sklar® US Army Retractor	60-1085
1	Sklar® Richardson Retractor	60-1676

41 Total pieces in this set (25 distinct instruments)

QTY	INSTRUMENT	SKLAR REF.
1	Sklar® Richardson Retractor	60-1681
1	Sklar® Zalkind Ribbon Retractor - 7-3/4" x 1"	22-8075
2	Sklar® Senn Retractor - 6-1/4", Sharp Prongs	60-6162
1	Emesis Basin - 9-7/8" x 4-1/2" x 2-1/8", 26 fl oz	10-1699
1	Sklar® Rochester-Pean Hemostatic Forceps - 5-1/2", Curved	17-2255
1	E-Z LOCK INSTR STRING 6 X 2.5	10-1806
41	Total pieces in this set (25 distinct instruments)	



Sklar — A 130-year heritage of surgical excellence

Sklar Corporation has been crafting surgical instruments in West Chester, Pennsylvania since 1882. Every instrument is forged from high-grade German stainless steel, finished by hand, and verified by Sklar's documented quality system before shipment.



51-point incoming inspection

Each component is checked against Sklar's quality standards before being approved for assembly.



6-point final inspection

Every assembled instrument is hand-checked against a final 6-point checklist before being shipped.

DE

German stainless steel

Surgical-grade material on OR Quality and Sklarlite™ lines. Hand-finished to exacting tolerances.



130+ years heritage

Crafting surgical instruments since 1882. Trusted by U.S. hospitals, surgery centers, and clinics.

Set reference data

SKLAR REFERENCE	98-1331
SET TYPE	Minor Basic Set — pre-configured (41 pieces)
QUALITY TIER	Sklar OR Quality — German stainless
REUSABILITY	Reusable · non-sterile · steam autoclave compatible
INSPECTION	51-point incoming + 6-point final (per Sklar IFU)
WARRANTY · STAINLESS	Lifetime against manufacturing defects
WARRANTY · CARBIDE	5 years (sharpen / repair / replace)
MANUFACTURER	Sklar Corporation · West Chester, PA · est. 1882
DISTRIBUTION	PCB Medical — authorized Sklar distributor

WARRANTY & QUALITY

Backed by Sklar's Lifetime Guarantee

"All Sklar Instruments are guaranteed for life against manufacturing defects, provided that the instrument is used for its intended surgical purpose. Sklar will repair or replace defective instruments at no charge, at our discretion."

✓ Lifetime guarantee

✓ 5-yr carbide warranty

✓ 51-point QC

✓ Made by Sklar

