



SKLAR OR QUALITY · PRE-CONFIGURED SET

LAP BOWEL RESECTION SET

Pre-configured 11-piece instrument set by Sklar Corporation (ref. 98-1357) — German stainless steel, factory-inspected, backed by Sklar's Lifetime Guarantee against manufacturing defects. Distributed in your market by PCB Medical.

German Stainless Steel

Lifetime Guarantee

51-point QC

5-yr carbide warranty

Autoclave-ready

Set composition

Standard PCB Medical / Sklar 98-1357 configuration. Custom configurations available on request.

QTY	INSTRUMENT	SKLAR REF.
1	Skartech™ 5000 High Frequency Suction Control Valve - Single Trumpet Valve	31-2227
1	Skartech™ 5000 Insulated Electrode, With Suction - Spatula Electrode	31-2235
1	5000 TC 5MM CVD LFT NON-INS	31-2552
1	Skartech™ 5000 Laparoscopic Metzenbaum Scissors Complete Instrument - 13", Curved, Blunt/Blunt, 5mm Tip, with 'XC' Non-Ratcheting Handle	31-9010XC
1	Skartech™ 5000 Micro Scissors - Curved. Complete Instrument 33cm 5mm	31-9025XC
1	Skartech™ 5000 Dolphin Forceps - w/ Spoon. Complete Instrument 33cm 5mm	31-9040XC
1	Skartech™ 5000 Grasping Forceps - Tapered, Rounded Jaws, Ratchet Handle, Complete Instrument 33cm 5mm	31-9051YC
1	Skartech™ 5000 Laparoscopic Allis TC Forceps Complete Instrument - 13", 5mm Tip, Atraumatic Cross-Serrated, with 'YC' Ratcheting Handle	31-9067YC
1	Skartech™ 5000 Grasping Forceps - w/ Spoon. Complete Instrument 33cm 5mm	31-9080YC
1	Skartech™ 5000 Fenestrated Grasping Forceps - Long Jaws. Complete Instrument 33cm 5mm	31-9095YC
1	Skartech™ 5000 Laparoscopic Babcock Forceps Complete Instrument - 13", 10mm Tip, with 'YC' Ratcheting Handle	31-9140YC
11	Total pieces in this set (11 distinct instruments)	



Sklar — A 130-year heritage of surgical excellence

Sklar Corporation has been crafting surgical instruments in West Chester, Pennsylvania since 1882. Every instrument is forged from high-grade German stainless steel, finished by hand, and verified by Sklar's documented quality system before shipment.



51-point incoming inspection

Each component is checked against Sklar's quality standards before being approved for assembly.



6-point final inspection

Every assembled instrument is hand-checked against a final 6-point checklist before being shipped.

DE

German stainless steel

Surgical-grade material on OR Quality and Sklarlite™ lines. Hand-finished to exacting tolerances.



130+ years heritage

Crafting surgical instruments since 1882. Trusted by U.S. hospitals, surgery centers, and clinics.

Set reference data

SKLAR REFERENCE	98-1357
SET TYPE	Lap Bowel Resection Set — pre-configured (11 pieces)
QUALITY TIER	Sklar OR Quality — German stainless
REUSABILITY	Reusable · non-sterile · steam autoclave compatible
INSPECTION	51-point incoming + 6-point final (per Sklar IFU)
WARRANTY · STAINLESS	Lifetime against manufacturing defects
WARRANTY · CARBIDE	5 years (sharpen / repair / replace)
MANUFACTURER	Sklar Corporation · West Chester, PA · est. 1882
DISTRIBUTION	PCB Medical — authorized Sklar distributor

WARRANTY & QUALITY

Backed by Sklar's Lifetime Guarantee

"All Sklar Instruments are guaranteed for life against manufacturing defects, provided that the instrument is used for its intended surgical purpose. Sklar will repair or replace defective instruments at no charge, at our discretion."

✓ Lifetime guarantee

✓ 5-yr carbide warranty

✓ 51-point QC

✓ Made by Sklar