



SKLAR OR QUALITY · PRE-CONFIGURED SET

# ABDOMINAL SET

Pre-configured 31-piece instrument set by Sklar Corporation (ref. 98-1400) — German stainless steel, factory-inspected, backed by Sklar's Lifetime Guarantee against manufacturing defects. Distributed in your market by PCB Medical.

German Stainless Steel

Lifetime Guarantee

51-point QC

5-yr carbide warranty

Autoclave-ready

## Set composition

Standard PCB Medical / Sklar 98-1400 configuration. Custom configurations available on request.

QTY	INSTRUMENT	SKLAR REF.
1	SklarLite™ Sterilization Container - Mid-Size, Standard Lid, Non-Perforated Bottom, 18-1/4" x 11" x 6", Silver	10-4394
1	SklarLite™ Sterilization Container Wire Basket - Mid-Size, 16" x 10" x 4"	10-4484
1	E-Z LOCK INSTR STRING 8 X 2.5	10-1808
1	Sklar Edge™ TC Mayo Dissecting Scissors - 9", Curved	16-1615
2	Sklar® Russian Tissue Forceps - 8"	19-2280
1	Sklar® Ochsner Trocar	34-1414
8	Sklar® Allis Tissue Forceps - 7-1/2", 5x6 Teeth	36-2275
1	Sklar® Jorgenson Scissors	52-2970
4	Sklar® Phaneuf Uterine Artery Forceps	90-2680
1	Sklar® O'Sullivan-O'Connor Abdominal Retractor	90-3025
2	Sklar® Sims Uterine Curette - 11", Sharp, #0	90-5500
2	Sklar® Barrett Tenaculum Forceps - 7"	91-1070
2	Sklar® Schroeder Tenaculum Forceps - 10-1/4", Round-Pointed Tips (Pozzi Style)	91-1298
2	Sklar® Staude-Moore Tenaculum Forceps	91-1465
2	Sklar® Jacobs Vulsellum Forceps	91-2282
<b>31</b>	<b>Total pieces in this set (15 distinct instruments)</b>	



## Sklar — A 130-year heritage of surgical excellence

Sklar Corporation has been crafting surgical instruments in West Chester, Pennsylvania since 1882. Every instrument is forged from high-grade German stainless steel, finished by hand, and verified by Sklar's documented quality system before shipment.



### 51-point incoming inspection

Each component is checked against Sklar's quality standards before being approved for assembly.



### 6-point final inspection

Every assembled instrument is hand-checked against a final 6-point checklist before being shipped.

### DE

#### German stainless steel

Surgical-grade material on OR Quality and Sklarlite™ lines. Hand-finished to exacting tolerances.



### 130+ years heritage

Crafting surgical instruments since 1882. Trusted by U.S. hospitals, surgery centers, and clinics.

## Set reference data

SKLAR REFERENCE	98-1400
SET TYPE	Abdominal Set — pre-configured (31 pieces)
QUALITY TIER	Sklar OR Quality — German stainless
REUSABILITY	Reusable · non-sterile · steam autoclave compatible
INSPECTION	51-point incoming + 6-point final (per Sklar IFU)
WARRANTY · STAINLESS	Lifetime against manufacturing defects
WARRANTY · CARBIDE	5 years (sharpen / repair / replace)
MANUFACTURER	Sklar Corporation · West Chester, PA · est. 1882
DISTRIBUTION	PCB Medical — authorized Sklar distributor

### WARRANTY & QUALITY

## Backed by Sklar's Lifetime Guarantee

*"All Sklar Instruments are guaranteed for life against manufacturing defects, provided that the instrument is used for its intended surgical purpose. Sklar will repair or replace defective instruments at no charge, at our discretion."*

✓ Lifetime guarantee

✓ 5-yr carbide warranty

✓ 51-point QC

✓ Made by Sklar

