



SKLAR OR QUALITY · PRE-CONFIGURED SET

# MAJOR PLASTIC SET

Pre-configured 127-piece instrument set by Sklar Corporation (ref. 98-1453) — German stainless steel, factory-inspected, backed by Sklar's Lifetime Guarantee against manufacturing defects. Distributed in your market by PCB Medical.

German Stainless Steel

Lifetime Guarantee

51-point QC

5-yr carbide warranty

Autoclave-ready

## Set composition

Standard PCB Medical / Sklar 98-1453 configuration. Custom configurations available on request.

QTY	INSTRUMENT	SKLAR REF.
2	Sklar® Scalpel Handle - #3	06-2903
1	Sklar® Scalpel Handle - #7	06-2907
1	Sklar® Weinstein Instrument Holder	10-2010
1	Skларcut™ Mayo Dissecting Scissors - 6-3/4", Curved, 1 Standard Edge, 1 Micro-Serrated Edge	15-3565
1	Skларcut™ Metzenbaum-Lahey Dissecting Scissors - 5-3/4", Curved, Blunt/Blunt, 1 Standard Edge, 1 Micro-Serrated Edge	15-3570
1	Skларcut™ Iris Scissors - 4-1/2", Straight, Sharp/Sharp, 1 Standard Edge, 1 Micro-Serrated Edge	15-3585
1	Skларcut™ Iris Scissors - 4-1/2", Curved, Sharp/Sharp, 1 Standard Edge, 1 Micro-Serrated Edge	15-3590
1	Skлар Edge™ TC Mayo Dissecting Scissors - 6-3/4", Straight	16-1510
1	Skлар Edge™ TC Metzenbaum-Lahey Dissecting Scissors - 5-3/4", Curved, Blunt/Blunt, Delicate	16-1905
1	Skлар Edge™ TC Metzenbaum Dissecting Scissors - 7", Curved, Blunt/Blunt, Delicate	16-2105
20	Sklar® Micro Mosquito Forceps	17-1225
4	Sklar® Halsted Mosquito Hemostatic Forceps - 5", Straight	17-1450
4	Sklar® Rochester-Pean Hemostatic Forceps - 6-1/4", Curved	17-2262
4	Sklar® Rochester-Ochsner Forceps - 6-1/4", Straight, Serrated, 1x2 Teeth	17-2862
4	Sklar® Crile Hemostatic Forceps - 5-1/2", Straight	17-3055
20	Sklar® Crile Hemostatic Forceps - 5-1/2", Curved	17-3155
<b>127</b>	<b>Total pieces in this set (46 distinct instruments)</b>	

QTY	INSTRUMENT	SKLAR REF.
2	Sklar® Dressing Forceps - 6", Straight, Serrated	19-1060
2	Sklar® Tissue Forceps - 6", Straight, 1×2 Teeth	19-1260
2	Sklar® TC Halsey / Webster Needle Holder - 5", Smooth	21-8000
2	Sklar® TC Webster / Halsey Needle Holder - 5-1/4", Smooth	21-8002
2	Sklar® TC Mayo-Hegar Needle Holder - 6", Serrated	21-8031
2	Sklar® TC Mayo-Hegar Needle Holder - 7", Serrated	21-8032
1	Sklar® Weitlaner Retractor - Sharp 3×4 Teeth 5-1/2"	22-8255
2	Sklarlite XD™ Mixer Right Angle Forceps	23-2521
1	Sklar® Director and Tongue Tie - 6"	30-1060
1	Sklar® Probe - 6", with Eye	30-1960
4	Sklar® Allis Tissue Forceps - 6", 5×6 Teeth"	36-2260
2	Sklar® Adson-Brown Tissue Forceps - 4-3/4", Straight, Side-Grasping 7×7 Teeth	47-1447
1	Sklar® Freer Elevator	47-2382
2	Sklar® Joseph Skin Hook - 6-1/4", 2 Prongs, Sharp, 2mm	47-2502
2	Sklar® Gillies Hook Large 4mm 7"	47-2572
2	Sklar® Blair-Rollet Retractor - 5-1/2", 4 Prongs, Sharp, Delicate	47-2655
4	Sklar® Lorna (Edna) Towel Clamp - 5-1/4"	47-2852
2	Sklar® Cushing Nerve and Vein Retractor	50-1290
2	Sklar® Adson Hemostatic Forceps - 7-1/4", Curved, Half-Serrated Jaws	50-1472
2	Sklar® Adson Tissue Forceps - 4-3/4", Straight, 1×2 Teeth	50-3047
2	Sklar® TC Adson Dressing Forceps - 4-3/4", Straight, Serrated	51-3447
2	Sklar® DeBakey Atraumatic Tissue Forceps	52-5177
2	Sklar® Richardson Retractor	60-1676
2	Sklar® Richardson Retractor	60-1681
2	Sklar® Volkman Retractor (USA) - 8-1/2", Sharp, 4 Prongs	60-5604
2	Sklar® Senn Retractor - 6-1/4", Sharp Prongs	60-6162
1	Sklar® Stevens Tenotomy Scissors - 4-1/2", Straight, Blunt/Blunt	64-1242
1	Sklar® Stevens Tenotomy Scissors - 4-1/2", Curved, Blunt/Blunt	64-1342
2	Sklar® Desmarres Lid Retractor - 6-1/4", 14mm	65-1842
4	Sklar® Foerster Sponge Forceps - 9-1/2", Straight, Serrated	87-2195

**127 Total pieces in this set (46 distinct instruments)**



## Sklar — A 130-year heritage of surgical excellence

Sklar Corporation has been crafting surgical instruments in West Chester, Pennsylvania since 1882. Every instrument is forged from high-grade German stainless steel, finished by hand, and verified by Sklar's documented quality system before shipment.



### 51-point incoming inspection

Each component is checked against Sklar's quality standards before being approved for assembly.



### 6-point final inspection

Every assembled instrument is hand-checked against a final 6-point checklist before being shipped.

### DE

#### German stainless steel

Surgical-grade material on OR Quality and Sklarlite™ lines. Hand-finished to exacting tolerances.



### 130+ years heritage

Crafting surgical instruments since 1882. Trusted by U.S. hospitals, surgery centers, and clinics.

## Set reference data

SKLAR REFERENCE	98-1453
SET TYPE	Major Plastic Set — pre-configured (127 pieces)
QUALITY TIER	Sklar OR Quality — German stainless
REUSABILITY	Reusable · non-sterile · steam autoclave compatible
INSPECTION	51-point incoming + 6-point final (per Sklar IFU)
WARRANTY · STAINLESS	Lifetime against manufacturing defects
WARRANTY · CARBIDE	5 years (sharpen / repair / replace)
MANUFACTURER	Sklar Corporation · West Chester, PA · est. 1882
DISTRIBUTION	PCB Medical — authorized Sklar distributor

### WARRANTY & QUALITY

## Backed by Sklar's Lifetime Guarantee

*"All Sklar Instruments are guaranteed for life against manufacturing defects, provided that the instrument is used for its intended surgical purpose. Sklar will repair or replace defective instruments at no charge, at our discretion."*

✓ Lifetime guarantee

✓ 5-yr carbide warranty

✓ 51-point QC

✓ Made by Sklar

