



SKLAR OR QUALITY · PRE-CONFIGURED SET

TRACHEOTOMY SET

Pre-configured 32-piece instrument set by Sklar Corporation (ref. 98-1468) — German stainless steel, factory-inspected, backed by Sklar's Lifetime Guarantee against manufacturing defects. Distributed in your market by PCB Medical.

German Stainless Steel

Lifetime Guarantee

51-point QC

5-yr carbide warranty

Autoclave-ready

Set composition

Standard PCB Medical / Sklar 98-1468 configuration. Custom configurations available on request.

QTY	INSTRUMENT	SKLAR REF.
1	Sklar® Scalpel Handle - #3	06-2903
1	Sklar® Yankauer Suction Tube - 12-1/4"	07-1820
1	Sklar® Mayo Dissecting Scissors - 5-1/2", Straight	15-1555
3	Sklar® Halsted Mosquito Hemostatic Forceps - 5", Straight	17-1450
3	Sklar® Halsted Mosquito Hemostatic Forceps - 5", Curved	17-1550
2	Sklar® Rochester-Pean Hemostatic Forceps - 5-1/2", Curved	17-2255
4	Sklar® Crile Hemostatic Forceps - 5-1/2", Curved	17-3155
1	Sklar® TC Crile-Wood Needle Holder - 6", Serrated, Left-Handed	21-9030
1	Sklar® Spencer (Littauer) Suture Scissors - 4-1/2"	22-2845
1	Sklar® Weitlaner Retractor - Sharp 3×4 Teeth 5-1/2"	22-8255
2	Sklar® Allis Tissue Forceps - 6", 5×6 Teeth"	36-2260
1	Sklar® Metzenbaum-Lahey Dissecting Scissors - 5-3/4", Curved, Blunt/Blunt	47-1360
4	Sklar® Backhaus Towel Clamp - 3-1/2"	47-2935
2	Sklar® Cushing Nerve Hook - #1, 7-1/2"	50-1307
1	Sklar® Adson Tissue Forceps - 4-3/4", Straight, 1×2 Teeth	50-3047
1	Sklar® Adson (Ewald) Dressing Forceps - 4-3/4", Straight, Serrated	50-3147
1	Sklar® Ballenger Sponge Forceps - 7", Curved, Serrated	74-1670
1	Sklar® News Tracheal Hook	78-3545
1	Sklar® Lukens Thymus Retractor	78-3580

32 Total pieces in this set (19 distinct instruments)



Sklar — A 130-year heritage of surgical excellence

Sklar Corporation has been crafting surgical instruments in West Chester, Pennsylvania since 1882. Every instrument is forged from high-grade German stainless steel, finished by hand, and verified by Sklar's documented quality system before shipment.



51-point incoming inspection

Each component is checked against Sklar's quality standards before being approved for assembly.



6-point final inspection

Every assembled instrument is hand-checked against a final 6-point checklist before being shipped.

DE

German stainless steel

Surgical-grade material on OR Quality and Sklarlite™ lines. Hand-finished to exacting tolerances.



130+ years heritage

Crafting surgical instruments since 1882. Trusted by U.S. hospitals, surgery centers, and clinics.

Set reference data

SKLAR REFERENCE	98-1468
SET TYPE	Tracheotomy Set — pre-configured (32 pieces)
QUALITY TIER	Sklar OR Quality — German stainless
REUSABILITY	Reusable · non-sterile · steam autoclave compatible
INSPECTION	51-point incoming + 6-point final (per Sklar IFU)
WARRANTY · STAINLESS	Lifetime against manufacturing defects
WARRANTY · CARBIDE	5 years (sharpen / repair / replace)
MANUFACTURER	Sklar Corporation · West Chester, PA · est. 1882
DISTRIBUTION	PCB Medical — authorized Sklar distributor

WARRANTY & QUALITY

Backed by Sklar's Lifetime Guarantee

"All Sklar Instruments are guaranteed for life against manufacturing defects, provided that the instrument is used for its intended surgical purpose. Sklar will repair or replace defective instruments at no charge, at our discretion."

✓ Lifetime guarantee

✓ 5-yr carbide warranty

✓ 51-point QC

✓ Made by Sklar

